



Utilizing nano-meter resolution laser Doppler vibrometer to measure small vibration of human skin

運用奈米精度雷射都卜勒測振儀偵測皮膚震動位移

1

Speaker: Yi-Cheng Tsai

Advisor: How-Foo Chen



2

Laser Doppler Technology Application Fields

Civil Engineering



Mechanical



Aerospace



Automotive



Biomedical





Laser Doppler Vibrometer

Benefit


- Remote (non-contact)
- High spatial resolution
- Reduced testing time
- Wide frequency range (near-DC to 1GHz)
- Large velocity range ($\sim 30\text{m/s}$)





Laser Doppler Principle



Doppler Effect



Doppler frequency shift


$$\Delta f = f'' - f_{laser} = \frac{\pm 2v_{sample}}{\lambda_{laser}}$$

Heterodyne Detection



Acousto-Optical
Modulator
(AOM)

Data Analysis





Hilbert-Huang Transform (HHT)

non-stationary
nonlinear data

Empirical Mode
Decomposition (EMD)

Intrinsic Mode
Function (IMF)


Experimental Setup





Experimental Design

Reference



Sample



Human skin surface physiological signals

Pulse Test Pulse at the wrist

Newborn : 130 - 150 bpm

Children : 65 - 105 bpm

Adults : 60 - 100 bpm

Test Object 25-year-old adult


Estimated pulse rate range

Adults : 60 - 100 bpm
1 - 1.67 Hz




7

Data Analysis - Hilbert-Huang Transform (HHT)



Matlab



Hilbert Transform
Instantaneous Frequency

Sample Surface
Displacement

Empirical Mode
Decomposition (EMD)

Intrinsic Mode
Function (IMF)

Find results
what we expected



Hilbert Spectral Analysis



Data Analysis

After Empirical Mode Decomposition (EMD)



Each Intrinsic Mode Function (IMF)





Data Analysis

Frequency Domain



Pulse
1.2774 Hz → 77 bpm


Adults : 60 - 100 bpm
1 - 1.67 Hz



10

Discussion

Calculate the Ratio




Calculate the ratio of the area of the main peak to the entire spectrum.

Ratio = 0.3284




Discussion


Pulse test




$R=0.3343$
 $F=1.0014$




$R=0.2314$
 $F=1.4305$




$R=0.2645$
 $F=1.1444$



$R=0.2445$
 $F=1.1444$



$R=0.2121$
 $F=1.0490$



$R=0.3603$
 $F=1.0967$




Measuring
position




Discussion


Not on pulse




R=0.1012
F=1.2398




R=0.1052
F=1.5259




R=0.0801
F=1.6689




R=0.1007
F=1.7166



R=0.0896
F=1.9074



R=0.0871
F=1.8120




Measuring
position




Discussion

Comparison


On pulse




$R=0.3343$
 $F=1.0014$




$R=0.2645$
 $F=1.1444$




$R=0.2121$
 $F=1.0490$



$R=0.2314$
 $F=1.4305$




$R=0.2445$
 $F=1.1444$




$R=0.3603$
 $F=1.0967$


Not on pulse




$R=0.1012$
 $F=1.2398$




$R=0.1052$
 $F=1.5259$




$R=0.0801$
 $F=1.6689$



$R=0.1007$
 $F=1.7166$



$R=0.0896$
 $F=1.9074$



$R=0.0871$
 $F=1.8120$



Discussion

Comparison of the pulse and non-pulse signals

- On the pulse, the ratio can reach 0.2121 to 0.3603.
- On non-pulse, the ratio is only 0.0801 to 0.1052.

In these two cases, the ratio differs by 2-4 times.

Feasibility

- From the results of data analysis, it verifies the feasibility of physiological signal measurement.



Thank You for listening